## EN A. V. Fistula Sets for Hemodialysis

For safety reasons, please read the instructions for use before using this device.

After reading the instructions, keep this manual in a convenient place for future reference.

This device is sterilized by gamma-ray irradiation.

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
### EN

- ① Needle protector
- ② Needle
- ③ Wing (with Needle hub)\*
- ④ Tubing
- ⑤ Clamp
- ⑥ Luer connector
- ⑦ Cap
- ⑧ Anti-needle stick device

\*Needle hub is component for the device with rotatable wing only.

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### EN

#### Intended use

This device is a disposable medical device, intended to be used as vascular access for hemodialysis.

#### Indication

This device is designed as a disposable medical product for hemodialysis patients.

#### Contraindication

Non-reusable

#### Warning

1. This device is designed for single use. Do not re-sterilize. Re-use of single-use devices creates a potential risk of biological contamination and/or impairment of functional capability. This may lead to adverse reactions, including severe illness or death of the patient or user.

#### Cautions

1. This device shall be used by a skilled physician or medical professional under the supervision, instructions and guidance of a skilled physician.
2. The tubing may kink during shipment or installation. Inspect the tubing for kinks prior to and during patient use. Do not use severely kinked tubing as hemolysis may occur.
3. Aseptic technique is required during use of the device.
4. Open the pouch to take out this device immediately before use.
5. Use this device only with a blood circuit which has locking port connectors.
6. Ensure all connections are secure and fluid-tight before use, and monitor regularly for leaks during patient use.
7. Hemolysis may occur at elevated negative or positive pressures. Check access and for kinks whenever excessive pressure is noticed. Stenosis of the access and improper placement of access needles will cause reduced blood flow which results in excess pressure.
8. Do not alter the device.
9. SB-KAWASUMI LABORATORIES, INC. assumes no responsibility for use of any product code with an incompatible blood circuit.
10. This device may crack when exposed to some medicinal agents. Be very careful not to allow the device to come into contact with disinfectants like alcohol, local anesthetics and fat emulsifying agents.
11. If any serious incident occurred in relation to this device, the user should report it to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.
12. Consult your physician immediately whenever an adverse reaction occurs.
13. The package must be visually inspected before use. Do not use if the sterile package is damaged or unintentionally opened before use.
14. This device contains di (2-ethylhexyl) phthalate (DEHP). Animal experiments have shown DEHP to be potentially toxic to reproduction. Given the current state of scientific knowledge, a risk (in case of long-term use) especially for premature male babies cannot be completely excluded. The application of medical products containing DEHP should be restricted to short term use for pregnant women, breastfeeding mothers, infants, and children.

#### Storage


To store this device, avoid direct sunlight, ultraviolet rays, high temperature, high humidity, water leaks, dust and severe vibration. This device must be stored in a clean, dry place at room temperature.

#### Discarding

After use, this device should be treated as medical waste. When using the device without an Anti-needle stick device, pay attention to the tip of the needle when discarding this device. Dispose of all used materials into an appropriate sharps container or according to facility protocol.


#### Directions for Use

1. Clean and disinfect the puncture site while taking care to avoid contaminating it.
2. When the device is removed from the package, it should be inspected for any damage.
3. Re-tighten the cap. And if this device has a clamp, close the clamp.
4. Remove the needle protector while avoiding contact contamination of the needle.
5. Fold the wings together securely with the bevel facing up. (Refer to the physician's prescription for bevel direction.)  
For rotatable wing type, dot marks are indicated on the needle hub to indicate bevel position. A black dot on the needle hub indicates the bevel is facing up and a red dot indicates the bevel is facing down. Check the position of the needle with these dot marks.
6. Insert the needle smoothly into the blood vessel, and gently guide it to achieve secure placement.
7. Observing blood pulsation, fill the tubing with blood.
8. Secure the needle alignment and fix the wing with tape.
9. Detach the cap and connect with a blood circuit aseptically.
10. To make a secure patient connector/fistula needle connection (See Figure. 1):  
a. Establish male/female luer connections first. Before engaging the locking cap, holding only the clear end, firmly push the male patient connector into the female fistula connector with a slight twisting motion. This connection must be secured before engaging the locking cap.  
b. Engage the red or blue locking cap for additional security.  
c. Check the connections regularly.
11. Needle positions and all connections must be checked carefully prior to dialysis and during dialysis. Connections and puncture sites should be visually inspected to detect leaks and avoid blood loss.

#### Anti-needle stick device

Upon completion of a dialysis session, slide the Anti-needle stick device to the puncture site until the lower jaw is positioned under the wings.

Under wings

Place the middle finger on styptic cotton and hook the index finger on the Finger shield.

To remove the needle from the patient, grasp the tubing directly and pull it carefully until the wings are fastened by the gap and you hear a "click."

#### Warranty

SB-KAWASUMI LABORATORIES, INC. will replace any product that is proven to be defective. SB-KAWASUMI LABORATORIES, INC. does not accept liability for any loss or damage, direct, consequential or special, arising out of the application or use of or inability to use this device. Before using, the user shall determine the suitability of this device for the intended use and the user assumes all risk and liability whatsoever in connection therewith.

#### Description of symbols

	Manufacturer		Authorized representative in the European Union
	Medical Device		Importer in the European Union
	Sterilized using irradiation		Do not re-use
	Caution		Single sterile barrier system
	Product containing DEHP		Consult instructions for use
	Catalogue number		Batch code
	Date of manufacture Country of manufacture (Thailand)		Use-by date
	Fragile, handle with care		Keep dry
	Keep away from sunlight		Non-pyrogenic fluid path
	This way up		Do not use if package is damaged and consult instructions for use
	Do not re-sterilize		Unique Device Identifier
	Packaging unit		

WEB SITE ADDRESS: <https://www.sb-kawasumi.jp/e/business/IFU.html>

## SB-KAWASUMI LABORATORIES, INC.

3-25-4 Tonomachi, Kawasaki-ku, Kawasaki-shi,  
Kanagawa 210-8602, Japan  
Tel: +81-44-589-8070 Fax: +81-44-589-8910

CE 0123

**Medical Technology Promedt Consulting GmbH**

Ernst-Heckel-Straße 7, 66386 St. Ingbert, Germany

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